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Fraunhofer IPT. Whereas the CFRP needle is stiffer and stronger than metallic ones, it enables a completely artefactfree Magnetic Resonance Imaging. For the production of these needles the Micro-Pultrusion Process was applied and brought to serial application. At an outside diameter of merely 1 200 μm , the needle offers three lumen with 272 μm , 312 μm and 612 μm inside diameter, a cutting edge for the direct placement of the puncture needle into human tissue as well as the biocompatible surface coating. The three lumen enable the parallel use of operation instruments e.g. endoscope and laser combined to a rinsing functionality. Standardised ports like Luer-Lock adapters allow the coupling with commercially available medical technology components. Furthermore a puncture needle with one lumen was developed with an outer diameter of 800 μm and an inner lumen of 312 μm .

MULTIFUNCTIONAL CFRP PUNCTURE NEEDLE FOR MRI-INTERVENTIONS

Your Benefit

For gentle minimally invasive diagnostics and therapy Magnet Resonance Imaging (MRI) enables high resolution exposures of soft tissue combined with a low radiation level. Thereby, metallic materials disturb the magnetic field which inhibits an exact imaging and terminates a successful minimally invasive therapy. Compared to metallic materials, the use of Carbon Fiber Reinforced Plastics (CFRP) enables mechanical properties such as strength and stiffness and guarantees an artefactfree imaging.

The advantage in using CFRP for minimally invasive surgery instruments is the highest possible accuracy of positioning the device relative to vulnerable structures in soft tissue. Therefore the variable characteristics and high mechanical properties of fiber reinforced plastics allow new design aspects for diverse kinds of medical devices.

Our Solution

A multifunctional puncture needle made of fiber reinforced plastics was developed at

Our Services

- Development and validation of the production technology and transfer to serial application of the Micro-Pultrusion Process
- Build-up and development of manufacturing systems for the manufacturing of Fiber Reinforced Plastic parts for medical applications
- Design and development of prototypes
- Execution of all process relevant and regulatory adaptations and qualifications
- Supervision of any kind of mechanical, biological and medical material analysis